




# BR A-FIB

## PATIENT INSTRUCTIONS



### SYMPTOM RECORDING

1. When you feel a symptom press either **RECORD / SEND** button.
2. One beep indicates monitor is recording.
3. To **SEND** the **EKG, REMOVE THE PATIENT LEAD** from monitor and call 800.590.5559 or (516.872.8001). Tell the technician your name.
4. Put phone mouthpiece over hole in front of monitor when instructed to by technician. Then press the **RECORD / SEND** button on side of monitor. When monitor stops making noise, pick up and speak with technician.
5. To erase the monitor memory, wait until the technician tells you it's OK, then insert the patient lead back into the monitor.

### SEND EKG TRANSMISSION TO:

1.800.590.5559

In Nassau 516.872.8001

If You Have Any Questions, Please Call 516.872.7001

### DIRECTIONS FOR RETURNING EVENT/LOOP MONITORS

#### RETURNING TO DOCTORS

When you are finished with the monitor, please take it and any leftover electrodes, batteries and equipment back to your doctor's office. The monitor can be returned to any nurse. You do not need to have or make an appointment for this. In the event that you are unable to take it back to your doctor, please contact the monitor return department at 516.872.7001. We will provide you with a UPS return envelope.

#### RETURNING TO SKY MEDICAL

When you are finished with the monitor, please place it and any leftover electrodes, batteries, and equipment into the UPS return envelope. Seal up the open end of the envelope and place it into any UPS Drop Off Box. You will find a list of Drop Off locations with the return envelope.

**THANK YOU**

