




# GEN PATIENT INSTRUCTIONS

## 3-WIRE ELECTRODE PLACEMENT



### SYMPTOM RECORDING

1. When you feel a symptom, push the red **STORE** button. A confirmation beep should be heard.
2. As soon as possible, call Toll Free 1.800.590.5559. (In Nassau call 516.872.8001) Tell the technician your name.
3. To send data, place the telephone mouthpiece over the recorder and then press the green **SEND** button.  
**Leave connected to telephone until sound stops.  
Pick up phone and talk to Technician.**
4. To **ERASE** memory, open battery compartment on left side for 10 seconds.

### SEND EKG TRANSMISSION TO:

1.800.590.5559

In Nassau 516.872.8001

If You Have Any Questions, Please Call 516.872.7001

## DIRECTIONS FOR RETURNING EVENT/LOOP MONITORS

### RETURNING TO DOCTORS

When you are finished with the monitor, please take it and any leftover electrodes, batteries and equipment back to your doctor's office. The monitor can be returned to any nurse. You do not need to have or make an appointment for this. In the event that you are unable to take it back to your doctor, please contact the monitor return department at 516.872.7001. We will provide you with a UPS return envelope.

### RETURNING TO SKY MEDICAL

When you are finished with the monitor, please place it and any leftover electrodes, batteries, and equipment into the UPS return envelope. Seal up the open end of the envelope and place it into any UPS Drop Off Box. You will find a list of Drop Off locations with the return envelope.

**THANK YOU**

